

I-MOVE

Patient transport chair

DEDICATED SOLUTIONS FOR YOUR PATIENTS' TRANSFERT & WELLBEING





Patient transport chair I-MOVE
For the patients' reception



Stackable patient transport chair I-MOVE
for optimized tidying

I-MOVE

DESCRIPTION

I-MOVE is a transport chair for your patients' reception and sitting transport.

Manageable, robust and attractive, it is now stackable. It is the modern and colored object for your department.

With our big range of accessories, you can adapt your use.

Overall width: 75 cm.

Overall length: 105 cm.

Overall height: 100 cm.

Weight: 35 kg.

2 years warranty.

Technical data I-Move

Steel frame epoxy coated.

Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person.

Front swivel wheels Ø 125 mm.

Cushion in foam width 5 cm.

Seamless finish, vinyl coating.

Color range: fuchsia, green, blue.

Dimensions seat and back: 46 x 46 cm

Seat height : 55 cm.


Sliding and self-blocking foot rest.

Polyurethan armrests, can be folded up.

Safe Working load : 200 kg.



Two positions leg rest (optional)



Different colors



Sliding and self-blocking foot rest



I-MOVE



Z.I. de Fétan - CS 11001
01604 TREVoux Cedex
FRANCE

Tél : +33 (0)4 74 08 71 71
Fax : +33 (0)4 74 08 71 72
E-mail : info@bmb-medical.com
www.bmb-medical.com

© 2017 -  - Tél. 06 77 96 82 22 / www.krea-communication.fr

